

AntiStick®

Radial fiber with AntiStick technology

Precise without exception

- Finest and most flexible fiber on the market
- Optimal solution for short segments
- Unique heat radiation for less sticking

With its thickness of just 1.0 mm, the AntiStick® fiber is the finest and most flexible on the market for the ablation of varicose veins.

The AntiStick® fiber in combination with the **simLa® 6** at the wavelength of 1940 nm enables a fluent and precise working due to a low tendency to sticking.



Ordering information

	REF	Fiber length	Introducer	Shipping unit
AntiStick® radial fiber (ASR) 400µm	270127S	260 cm	4F / 17G	5 pcs
AntiStick® radial fiber (ASR) 600µm	270128S	260 cm	6F	5 pcs

iMS GmbH
innovative Medical Solutions
Riedstr. 63 A
82327 Tutzing
Germany
www.ims-medical.de

simLa® and AntiStick® are
registered trade marks
simLa® is protected by patent

Product listed may not be available
in all markets.